

# DRE Milano OB20

## OB/GYN Procedure Chair

### Electrically operated chair with modern European style

The New DRE Milano OB20 Procedure Chair features full electric operation. It has a European design and the surface is thermoplastic over a steel frame for easy cleaning. The Milano OB20 provides simple, automatic operation with 3 independent motors. Perfectly suited for use in gynecology, urology, colposcope exams, diagnostics and ultrasounds.

#### FEATURES

- Automatic reset switch that returns the table to its starting position
- 20" minimum height ensures easy access for the patient
- 40" maximum height ensures optimum position for the specialist
- Steel frame covered with thermoplastic for a sanitary and attractive appearance
- Seamless anatomical mattress for easy sanitation
- 3 independent electric movements: height, back rest and Trendelenburg
- 90° rotation to both sides
- Ergonomic styling
- Adjustable headrest



#### SPECIFICATIONS



**Maximum Patient Weight:**  
485 lbs. (200kg)



**Chair Length:**  
51 - 60in (129 - 151cm)

**Chair Width:**  
23in (59cm)

**Maximum Chair Height:**  
40in (102cm)

**Minimum Chair Height:**  
20in (52cm)



**Voltage:**  
230/110 V

**Frequency:**  
50/60 Hz

**Absorption:**  
1A max

**Transformer:**  
25V

**Phase:**  
L + N

#### STANDARD ACCESSORIES

- Removeable stainless steel bowl
- Roll holder
- Seat guides handles
- Fixed clamps
- Two thigh supports

381-2018-06-26