

An Elegant Solution for Today's Podiatry Practice

The DRE Milano P20 is a new take on traditional podiatry procedure chairs, outpacing other models in features, style and affordability. The Milano P20 comes standard with 180° chair rotation, unique split leg supports, hand remote and foot controls and standard armrests.

The Milano P20 is expertly designed for comfort for podiatrists and patients alike, comfortable for both routine visits and lengthy procedures. Using the Milano P20 podiatry chair's adjustable height, operators can lower the chair to 20.5 inches. The chair's many other physician-configured adjustments including the split leg controls, provide excellent visibility within the operative field. It includes a wide variety of power positions and a smart "auto return" feature. The chair will lay completely flat in a table position.

Features:

- 180° chair rotation: Unique swivel feature for convenient access to patient area.
- Leg supports swivel independently up to 50 degrees.
- Foot pedal adjusts the height, back and level positioning.
- 3 Electric movements: height, backrest and seat tilt.
- Partial Trendelenburg positioning.
- Foot supports adjusted with gas spring.
- Sleek system with a small footprint.
- Stitchless mattress/covering allows for easy cleanup and sterilization.



DRE Milano P20

Podiatry Chair

Equipment for the way *you* operate

Specifications:

- Voltage: 110V AC, 230V AC
- Frequency: 50/60 Hz
- Minimum chair height: 20" (52 cm)
- Maximum chair height: 32" (82 cm)
- Angle of Trendelenburg inclination: 30°
- Angle of backrest inclination: 80°
- Angle of leg rest inclination: 90°
- Mattress length: 66" (168 cm) | 73" (185 cm) with extended leg rest
- Mattress width: 23" (59 cm)
- Rotation: 90° on both sides

Optional Accessories:

- Foot supports extender cushion
- Hand set control
- Roll holder
- Foot tray
- Manicure bowl
- Alternate armrests
- Fabric leg rest covering
- Lamp adapter
- Wireless foot pedal

