

# DRE Milano T50

## Power Procedure Table

### Fully electric, multi-function power procedure table

The DRE Milano T50 is a fully electric, multi-function table that compares with popular brands in terms of style and functionality, at a fraction of the cost. With four actuators to adjust height, backrest, leg rest and Trendelenburg positioning, clinicians are ensured a clear path to the patient. The table's low height capabilities give the patient easy access in and out of the chair.

The padded, anatomically-configured mattress foam and adjustable or changeable headrest make the DRE Milano T50 Power Procedure Table naturally adaptable to a variety of lengthy office-based procedures. It also comes standard with hand and foot controls with smart return, allowing clinicians to easily maneuver the chair into a neutral starting position.

#### FEATURES

- Fully electric operation with four actuators to adjust table height, backrest, leg rest and Trendelenburg positioning.
- Standard handheld remote control and foot pedal include automatic reset function.
- Adjustable and removable headrest.
- Seat and back rails fit a variety of medical accessories.
- Magnetic back remote control easily attaches to the side of the chair.
- ABS coated steel frame gives the table a smooth finish while protecting the inner table components.
- Upholstery is designed for medical use and is resistant to water, air and disinfectants.
- Anatomic mattress is made of top-quality foam and includes a five-year warranty.
- Constructed using latex-free, non-toxic and non-carcinogenic materials.



*Shown with optional medical arm-rests*

#### SPECIFICATIONS



**Chair Weight:** 196 lbs. (89kg)

**Maximum Capacity Load:**  
440 lbs. (200kg)



**Adjustable Height:**  
24" (62cm) to 36" (92cm)

**Chair Dimensions:**  
74" x 24" x 34"  
(190 x 62 x 86cm)



**Trendelenburg:** 25°

**Backrest:** 75°

**Leg rest:** 65°



**Voltage:**  
110 and 220 V cc

#### COLOR OPTIONS



Blue



Cream



Chrome



Black

**WARRANTY:** 2 Years

390-2018-06-26

## Multi-Specialty Use

- › General Medical
- › Dermatology
- › Podiatry
- › Plastic, Cosmetic and Aesthetics, OB/GYN
- › ENT/Otolaryngology
- › Oral & Maxillofacial
- › Dental
- › Orthodontics

*Shown with optional accessories*

