

DRE Millennium 5

Hospital Stretcher

Available with easy-to-steer fifth wheel system

BENEFITS

- High-low adjustable mattress range is 13.5"
- Back-rest lifting range is 0 - 75° and controlled by oil hydraulics
- Knee-rest lifting range is 0 - 30° and controlled by adjustable crank
- Central braking system with steering function
- Trendelenburg range is $\pm 18^\circ$
- Two IV poles (one at the head and foot)
- 4" high density PU foam mattress covered with synthetic leather
- Collapsible side rails on both sides that can fold manually
- 8" steering wheel central brake on four casters

FEATURES

- 4" x 28.5" x 71.5" Pad
- 8" Casters
- 27" Wide Patient Surface
- Fold-Down Side Rails; Side Rails 14.5" x 39"
- Stationary Foot Section
- Four Brake/Steer Pedals
- Dual Pedestal Hydraulics
- Dual Side-Mounted Foot Control Hydraulics
- Hands-Free Trendelenburg and Reverse Trendelenburg
- Steel bed frame with powder coating
- Mattress base made of ABS plastic board
- Bumpers made of durable ABS plastic and located at each corner
- Base cover made of ABS plastic with space for an oxygen bottle and emergency devices



SPECIFICATIONS



Weight Capacity: 550 lbs



Overall Length: 83"

Overall Width: 33"

Mattress Surface:
71.5" (L) x 28.5" (W)

Height Adjustment Range:
22.5" - 36.5" with 8" Casters



Trendelenburg/
Reverse Trendelenburg: $\pm 18^\circ$

Backrest Adjustments: +90°



Casters:
8" with central locking one
tracking 4 brake/steer pedals

I.V. Attachment Sockets:
2

The fifth wheel
system saves
power and
allows for easy
maneuverability

397-2018-06-22