

DRE Waveline EZ MAX


Vital Signs Monitor

High-resolution color touch screen

The DRE Waveline EZ MAX is a full-featured touch-screen monitor designed for a variety of settings. The customizable menu allows operators to easily navigate between the different parameters. The DRE Waveline EZ MAX features a no-fan design, ideal for quiet recovery and sterile environments.

Features

- ✓ 10.4" touch-screen color TFT-LCD monitor
- ✓ 3/5-lead ECG, RESP, NIBP, SpO2, 2-TEMP, PR
- ✓ Arrhythmia, S-T Analysis, OxyCRG, Large Font Display
- ✓ USB data storage and review
- ✓ Bi-directional Communications with Central Monitoring System
- ✓ No-fan design
- ✓ Waveform color adjustable
- ✓ Customizable alarm setting



SPECIFICATIONS



Overall Dimensions:
12" L x 6.3" W x 10" H
(303 x 161 x 254mm)

Display:
10.4" touch screen color TFT,
800 x 600 resolution



Power:
100-240V, 50/60Hz



Parameters Measured:
SpO2, RESP,
NIBP, 2-TEMP,
3/5-lead ECG,
and PR

Optional Parameter

- ✓ Nellcor OxiMax™ SPO2
- ✓ 2-IBP
- ✓ ICP
- ✓ C.O. Respiration CO2 (Mainstream and Sidestream)

914-2017-11-14