

The PathAguard[®] ABLiS[®]. This FDA Class 1 Medical Device is a unique bathing system designed to mitigate the risk for Pathogens to grow.

A GLOBAL HEALTHCARE CHALLENGE

Healthcare associated infections acquired during hospital stays are the most common complications of hospital care and one of the most serious patient concerns.

The rise of antibiotic resistant pathogens coupled with increased incidence of Healthcare Associated Infections continues to put a significant strain on Healthcare infection control measures.

Introducing PathAguard[®] ABLiS[®]

The PathAguard[®] ABLiS[®] represents a welcome addition to your infection control arsenal. The PathAguard[®] ABLiS[®] is a unique patient bathing system designed to inhibit the spread of Pathogens that may cause HAI's. The reusable basin, is protected with silver ion technology which inhibits the growth of Pathogens providing product protection.


The single use liner is protected with silver ion technology providing an extra physical barrier, mitigating the opportunity for Pathogens to grow.

BENEFITS

- Inhibits the growth of Pathogens on the product
- Patient comfort
- High patient and healthcare provider satisfaction.
- Cost effective solution.

FEATURES

- Inhibits the growth of Pathogens on the products.
- FDA Class 1 Medical Device
- CE marked
- Latex-free single use liner
- Utilising Biomaster silver ion technology
- Tested and proven in a laboratory setting
- Ergonomically designed system
- Recyclable material



"In my estimation the ABLiS represents the cleanest, safest, and most effective approach to protecting against healthcare acquired infection during hygiene activity."


Dr Roy Sleator
BSc, MA, PhD, CBiol, FRSB
Microbiologist of the year 2016



The PathAguard[®] ABLiS[®] brings a new bed bathing experience to the healthcare setting. It is easy to use and offers a safe environment for both patient and healthcare professional.

The Biomaster silver ion technology works in a multi modal way inhibiting the ability of the Pathogens to grow on the products.

HOW DOES IT WORK?



EASY TO USE - 6 SIMPLE STEPS



1

Place liners inside Receptacle with liners feeding from top.



2

Pull liner until it reaches perforation and tear gently.



3

Press liner down to ensure snug fit.



4

Attach liner using prepunched holes over hooks in ABLiS.



5

Fill ABLiS to required level as per facility protocol



6

After use, dispose of contents and single use liner as per facility protocol.

INHIBITS THE GROWTH OF PATHOGENS ON THE PRODUCTS.

The PathAguard[®] ABLiS[®] has been tested for its efficacy in independent laboratory environments and consistently demonstrates that it inhibits the growth of Pathogens on the products.

Why use the ABLiS Liner with the basin?

"The design of the system by using a liner reduces patient exposure to cross-infection from reused basins and contributes to overall infection control strategy in any healthcare setting"

R. J. Russell

BA, Mod, PhD (Associate Professor of Microbiology)